

Product Datasheet

Rubella Virus E1 Mosaic Antibody (orb623312)

Description	Rubella Virus E1 Mosaic Antibody
Reactivity	Human
Tested Applications	ELISA, WB
Immunogen	Recombinant protein contains the Rubella Virus E1 regions, 157-176, 374-390, 213-239 amino acids
Form/Appearance	Liquid
Storage	Store at room temperature for up to 3 weeks. Store desiccated at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.