

Product Datasheet

Herpes Simplex Virus-1 gG Antibody (orb623315)

Description	Herpes Simplex Virus-1 gG Antibody
Reactivity	Human
Tested Applications	ELISA, WB
Immunogen	Recombinant protein contains the HSV-1 gG immunodominant regions, 84-175 amino acids.
Form/Appearance	Liquid
Storage	Store at room temperature for up to 3 weeks. Store desiccated at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.