

Product Datasheet

Human Angiostatin Kringles 1-3 protein (orb623348)

| | |
|------------------------|--|
| Description | Human Angiostatin Kringles 1-3 protein |
| Reactivity | Human |
| Form/Appearance | lyophilized |
| Storage | Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |