

## Product Datasheet

### Human BAI1-Associated 2 protein (orb623386)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Human BAI1-Associated 2 protein  |
| <b>Reactivity</b>      | Human  |
| <b>Form/Appearance</b> | Liquid   |
| <b>Storage</b>         | Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles. |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <a href="#">Q9UQB8</a>   |
| <b>Entrez</b>          | <a href="#">10458</a>  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |