

## Product Datasheet

### Human CD300C protein (orb623457)

<b>Description</b>	Human CD300C protein
<b>Reactivity</b>	Human
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage, it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q08708</a>
<b>Entrez</b>	<a href="#">10871</a>
<b>Expiration Date</b>	12 months from date of receipt.