

## Product Datasheet

### Endomorphin-2 Antibody (orb623488)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Endomorphin-2 Antibody   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse, Rat  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ICC, IF, IHC   |
| <b>Immunogen</b>           | Endomorphin-1 conjugated to BSA.   |
| <b>Form/Appearance</b>     | Whole serum  |
| <b>Storage</b>             | The product can be stored as supplied for up to 12 months at 2°C-4°C. After reconstitution, aliquot and store at -20°C for higher stability or at 4°C with an appropriate antibacterial agent. Avoid freeze thaw-cycles. |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | Rabbit IgG   |
| <b>Dilution Range</b>      | Immunohistochemistry: 5-10 ug/ml<br>Immunofluorescence: 5-15 ug/ml   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |