

Product Datasheet

LPAM-1 antibody (PE) (orb623588)

Description	Rat monoclonal antibody against LPAM-1 conjugated to PE
Reactivity	Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	TK1 cells
Target	LPAM-1
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Rat IgG2a kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.