
Product Datasheet

**SAE2/UBA2 Antibody (monoclonal, 5B13)
(orb623715)**

Description

SAE2/UBA2 Antibody (monoclonal, 5B13)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

E. coli-derived human SAE2/UBA2 recombinant protein (Position: E449-K564).

Form/Appearance

Lyophilized

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Add 0.2ml of distilled water will yield a concentration of 500µg/ml.

Isotype

Mouse IgG2b

Clonality


Monoclonal

MW


90 kDa

Uniprot ID
[Q9UBT2](#)
Expiration Date


12 months from date of receipt.



IHC analysis of SAE2/UBA2 using anti-SAE...



IHC analysis of SAE2/UBA2 using anti-SAE...



IHC analysis of SAE2/UBA2 using anti-SAE...