
Product Datasheet

Rhesus macaque Glypican 2 / GPC2 Protein (orb624213)

Description

Rhesus macaque Glypican 2, His Tag

(HEK293). It contains AA Ser 24...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Glypican 2

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 59.9 kDa. The protein migrates as 55-95 kDa under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

XP_001103286

Purity

90%

Source

Rhesus macaque Glypican 2, His Tag (orb624213) is expressed from human 293 cells (HEK293). It contains AA Ser 24 - Ser 552 (Accession # XP_001103286).

MW

59.9 kDa

NCBI[XP_001103286](#)**Expiration Date**

6 months from date of receipt.