


## Product Datasheet

### Human RBM24 ELISA kit (orb624248)

<b>Description</b>	Human RBM24 ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/ml
<b>Target</b>	RBM24
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9BX46</b>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
RBM24(RNA  
Binding  
Motif  
Protein  
24...