


Product Datasheet

Human LMNA ELISA kit (orb624270)

Description	Human LMNA ELISA kit
Reactivity	Human
Range	0.156-10ng/ml
Target	LMNA
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P02545
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
LMNA(Lamin-
A) ELISA Kit