


Product Datasheet

Human CLUSTERIN ANTIBODY ELISA kit (orb624289)

Description	Human CLUSTERIN ANTIBODY ELISA kit
Reactivity	Human
Range	1.563-100ng/ml
Target	CLUSTERIN ANTIBODY
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Human CLUSTERIN
ANTIBODY(CLUSTERIN
ANTIB...