


Product Datasheet

Human C5AR1 ELISA kit (orb624299)

| | |
|------------------------|--|
| Description | Human C5AR1 (C5a anaphylatoxin chemotactic receptor 1) ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/ml |
| Target | C5AR1 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Sample Types | Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples |
| Uniprot ID | P21730 |
| Sensitivity | 0.094ng/ml |
| Expiration Date | Please enquire. |



Human C5AR1(C5a anaphylatoxin chemotacti...