


## Product Datasheet

### Animal FSH ELISA kit (orb624320)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Animal FSH ELISA kit   |
| <b>Reactivity</b>      | Equine   |
| <b>Range</b>           | 3.12-200mIU/ml   |
| <b>Target</b>          | FSH  |
| <b>Storage</b>         | 2-8 °C for 6 months  |
| <b>Note</b>            | For research use only  |
| <b>Sample Types</b>    | Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples |
| <b>Uniprot ID</b>      | <b>P01226</b>  |
| <b>Sensitivity</b>     | 1.875mIU/ml  |
| <b>Expiration Date</b> | Please enquire.  |



Horse  
FSH(Follicle  
Stimulating  
Hormone) ...