


Product Datasheet

Mouse Dapk1 ELISA kit (orb624336)

Description	Mouse Dapk1 ELISA kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Dapk1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q80YE7
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Dapk1
(Death-
associated
protein ki...