


## Product Datasheet

### Porcine Triiodothyronine ELISA kit (orb624344)

<b>Description</b>	Porcine Triiodothyronine ELISA kit
<b>Reactivity</b>	Porcine
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	T3
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Porcine  
T3(Triiodothyronine)  
ELISA Kit