


Product Datasheet

Porcine Thyroxine ELISA kit (orb624345)

Description	Porcine Thyroxine ELISA kit
Reactivity	Porcine
Range	0.703-45ng/ml
Target	T4
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	0.422ng/ml
Expiration Date	Please enquire.



Porcine
T4(Thyroxine)
ELISA Kit