


Product Datasheet

Rat Oxytocin ELISA kit (orb624349)

Description	Rat Oxytocin ELISA kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	OT
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01179
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
OT(Oxytocin)
ELISA Kit