


Product Datasheet

Rat Substance P ELISA kit (orb624350)

Description	Rat Substance P ELISA kit
Reactivity	Rat
Range	7.813-500pg/ml
Target	SP
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P06767
Sensitivity	4.688pg/ml
Expiration Date	Please enquire.



Rat
SP(Substance
P) ELISA Kit