


## Product Datasheet

### 4-HNE ELISA kit (orb624366)

<b>Description</b>	4-HNE ELISA kit
<b>Reactivity</b>	All
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	4-HNE
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



4-HNE (4-Hydroxynonenal)  
ELISA Kit