

Product Datasheet

ITGB3 (phospho-Tyr773) antibody (orb6247)

Description	ITGB3 (phospho-Tyr773) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ITGB3 around the phosphorylation site of Tyr773 (PL(p-Y)KE)
Target	phospho-ITGB3 (Tyr773)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	87 kDa
Uniprot ID	P05106
Expiration Date	12 months from date of receipt.