

Product Datasheet

CACNA1S antibody (orb625347)

Description	Rabbit polyclonal antibody to CACNA1S
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	calcium channel, voltage-dependent, L type, alpha 1S subunit
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	200-220kd
Uniprot ID	Q13698
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.