



Product Datasheet

CTDSP2 antibody (orb626099)

Description	Rabbit polyclonal antibody to CTDSP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	CTD(carboxy-terminal domain, RNA polymerase II, polypeptide A) small phosphatase 2
Target	CTDSP2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	31 kDa
Uniprot ID	O14595
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



human kidney tissue were subjected to SD...



Immunohistochemistry of paraffin-embedde...