



Product Datasheet

EDAR antibody (orb626609)

Description	Rabbit polyclonal antibody to EDAR
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, IP, WB
Immunogen	ectodysplasin A receptor
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	48 kDa
Uniprot ID	Q9UNE0
Dilution Range	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 cells using...



Immunohistochemistry of paraffin-embedde...