

Product Datasheet

GFP tag antibody (orb627318)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to GFP tag |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, IP, WB |
| Immunogen | GFP tag |
| Form/Appearance | liquid |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.) |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥95% as determined by SDS-PAGE |
| MW | 30 kDa |
| Dilution Range | WB: 1:2000-1:16000, IP: 1:1000-1:4000, IF: 1:20-1:200 |
| Expiration Date | 12 months from date of receipt. |