

Product Datasheet

GST Tag antibody (orb627538)

Description	Rabbit polyclonal antibody to GST Tag
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant protein
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Dilution Range	WB: 1:2000-1:20000
Expiration Date	12 months from date of receipt.