

Product Datasheet

HA-Tag antibody (orb627612)


Description	Mouse monoclonal antibody to HA-Tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	HA tag
Target	HA-Tag
Preservatives	liquid
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Clone Number	6A3
Purity	≥95% as determined by SDS-PAGE
Dilution Range	WB: 1:5000-1:50000, IP: 1:2000-1:20000, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



HA-tagged fusion protein were subjected to SDS PAGE followed by western blot with orb627612 (HA tag Antibody) at dilution of 1:20000.



IP Result of anti-HA-tag (IP: orb627612, 5 ug; Detection: orb627612 1:5000) with Transfected HEK-293 cells lysate 400 ug.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)