

## Product Datasheet

### HA-Tag antibody (orb627612)

<b>Description</b>	Mouse monoclonal antibody to HA-Tag
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IP, WB
<b>Immunogen</b>	HA tag
<b>Form/Appearance</b>	liquid
<b>Storage</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥95% as determined by SDS-PAGE
<b>Dilution Range</b>	WB: 1:5000-1:50000, IP: 1:2000-1:20000, IF: 1:10-1:100
<b>Expiration Date</b>	12 months from date of receipt.