

Product Datasheet

IL-16 antibody (orb627993)

Description	Mouse monoclonal antibody to IL-16
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	interleukin 16(lymphocyte chemoattractant factor)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	37 kDa
Uniprot ID	Q14005
Dilution Range	WB: 1:5000-1:20000
Expiration Date	12 months from date of receipt.