

Product Datasheet

MBP tag antibody (orb628680)

Description	Mouse monoclonal antibody to MBP tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	MBP-Tag
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	40 kDa
Dilution Range	WB: 1:1000-1:8000, IP: 1:5000-1:50000
Expiration Date	12 months from date of receipt.