

Product Datasheet

LAT (phospho-Tyr171) antibody (orb6291)

Description	LAT (phospho-Tyr171) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human LAT isoform b around the phosphorylation site of Tyr171 (D(p-Y)VN)
Target	Phospho-LAT (Tyr171)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29 kDa
Uniprot ID	O43561
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.