

Product Datasheet

LAT (phospho-Tyr191) antibody (orb6292)

Description	LAT (phospho-Tyr191) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human LAT isoform b around the phosphorylation site of Tyr191 (RE(p-Y)VN)
Target	Phospho-LAT (Tyr191)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29 kDa
Uniprot ID	O43561
Expiration Date	12 months from date of receipt.