

## Product Datasheet

### NFX1 antibody (orb629284)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to NFX1  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, WB   |
| <b>Immunogen</b>           | nuclear transcription factor, X-box binding 1   |
| <b>Form/Appearance</b>     | liquid  |
| <b>Storage</b>             | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.) |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Purity</b>              | ≥95% as determined by SDS-PAGE  |
| <b>MW</b>                  | 124kd   |
| <b>Uniprot ID</b>          | <a href="https://www.uniprot.org/entry/Q12986">Q12986</a>   |
| <b>Dilution Range</b>      | WB: 1:500-1:2000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |