

Product Datasheet

Lck (phospho-Tyr505) antibody (orb6309)

Description	Lck (phospho-Tyr505) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Lck around the phosphorylation site of Tyr505 (GQ(P-Y)QP)
Target	Phospho-Lck (Tyr505)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	P06239
Expiration Date	12 months from date of receipt.