

Product Datasheet

UBC12 antibody (orb632291)

Description	Rabbit polyclonal antibody to UBC12
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	ubiquitin-conjugating enzyme E2M(UBC12 homolog, yeast)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60 kDa
Uniprot ID	P61081
Dilution Range	WB: 1:200-1:2000, IF: 1:20-1:200;IP: 1:200-1:2000
Expiration Date	12 months from date of receipt.