

Product Datasheet

VCPIP1 antibody (orb632497)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to VCPIP1 |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, WB |
| Immunogen | valosin containing protein(p97)/p47 complex interacting protein 1 |
| Form/Appearance | liquid |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.) |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% as determined by SDS-PAGE |
| MW | 89 kDa |
| Uniprot ID | Q96JH7 |
| Dilution Range | WB: 1:200-1:2000, IHC-P: 1:20-1:200 |
| Expiration Date | 12 months from date of receipt. |