



Product Datasheet

Cardiac Troponin I antibody (orb632856)


Description	Mouse monoclonal antibody to Cardiac Troponin I
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	troponin I type 3(cardiac)
Target	cTnI
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	24-26 kDa
Uniprot ID	P19429
Dilution Range	WB: 1:5000-1:50000, IHC-P: 1:50-1:500, IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human heart tissue.



Immunofluorescent analysis of mouse heart tissue.



Western blot analysis of Rat heart tissue.