

Product Datasheet

Streptavidin antibody (HRP) (orb633180)

Description	Rabbit polyclonal antibody against Streptavidin conjugated to HRP
Reactivity	Bacteria
Tested Applications	ELISA, WB
Preservatives	3mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Tag	N-terminal His Tag
Note	For research use only
Application notes	Centrifuge the vial prior to opening, reconstitute in sterile distilled water or buffer of interest to concentration of 0.1-1mg/ml by gently pipetting 2-3 times, don't vortex
Purity	Greater than 90% as determined by SDS-PAGE
Source	E.Coli
Uniprot ID	P22629
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:2000-1:10000 for western blotting with ECL substrates 1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)