
Product Datasheet

Affinipure Goat Anti-Rabbit IgG (H+L) antibody (HRP) (orb633183)

Description

Goat polyclonal antibody against

(H+L) conjugated to HRP

Species/Host

Goat

Reactivity

Rabbit

Conjugation

HRP

Tested Applications

ELISA, IHC, WB

Preservatives3mg/ml rAlbumin (IgG-Free,
Protease-Free)**Form/Appearance**

liquid

Concentration3mg/ml BSA (IgG-Free, Protease-
Free)**Storage**0.01M Sodium phosphate, 0.25M
NaCl, 50% glycerol, pH7.6. -20°C for
24 months , avoid repeated freezing
and thawing.**Note**

For research use only

Isotype

IgG

Clonality

Polyclonal

Hazard InformationThis product is for research use
only, not for diagnostic or
therapeutic use.**Dilution Range**1:1000-1:10000 for western blotting
with ECL substrates1:1000-1:20000
for ELISA and Western blotting with
chromogenic substrates1:20-1:50
for IHC**Expiration Date**

12 months from date of receipt.