

Product Datasheet

Donkey Anti- Goat IgG (H+L) antibody (HRP) (orb633184)

Description	Goat polyclonal antibody against Donkey Anti- Goat IgG (H+L) conjugated to HRP
Species/Host	Donkey
Reactivity	Goat
Conjugation	HRP
Tested Applications	ELISA, WB
Preservatives	3mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:2000-1:10000 for western blotting with ECL substrates 1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates
Expiration Date	12 months from date of receipt.