

Product Datasheet

Rabbit Anti-Sheep IgG (H+L) antibody (HRP) (orb633195)

Description	Rabbit polyclonal antibody against Rabbit Anti-Sheep IgG (H+L) conjugated to HRP
Species/Host	Rabbit
Reactivity	Sheep
Conjugation	HRP
Tested Applications	ELISA, WB
Preservatives	3mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Concentration	3mg/ml BSA (IgG-Free, Protease-Free)
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:2000-1:10000 for western blotting with ECL substrates 1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates
Expiration Date	12 months from date of receipt.