

Product Datasheet

Affinipure Donkey Anti- Rat IgG (H+L) antibody (AP) (orb633204)

Description	Rat monoclonal antibody against Affinipure Donkey Anti- Rat IgG (H+L) conjugated to AP
Species/Host	Donkey
Reactivity	Rat
Conjugation	AP
Tested Applications	ELISA, WB
Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Storage	0.01M Tris-HCl, 0.25M NaCl, 50% glycerol, pH8.0. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:1000-1:3000 for ELISA and Western blotting
Expiration Date	12 months from date of receipt.