

Product Datasheet

Affinipure Goat Anti Rabbit IgG (H+L) antibody (FITC) (orb633210)

Description	Rabbit polyclonal antibody against Affinipure Goat Anti Rabbit IgG (H+L) conjugated to FITC
Species/Host	Goat
Reactivity	Rabbit
Conjugation	FITC
Tested Applications	FC, IF
Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Concentration	5mg/ml BSA (IgG-Free, Protease-Free)
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol pH7.6. -20°C for 24 months.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:20-1:100 for most applications
Expiration Date	12 months from date of receipt.