

## Product Datasheet

# Affinipure Rabbit Anti-Goat IgG (H+L) antibody (FITC) (orb633212)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody against Affinipure Rabbit Anti-Goat IgG (H+L) conjugated to FITC |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Goat  |
| <b>Conjugation</b>         | FITC  |
| <b>Tested Applications</b> | FC, IF  |
| <b>Preservatives</b>       | 5mg/ml rAlbumin ( IgG-Free, Protease-Free)  |
| <b>Form/Appearance</b>     | liquid  |
| <b>Concentration</b>       | 5mg/ml BSA ( IgG-Free, Protease-Free)   |
| <b>Storage</b>             | 0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol pH7.6. -20°C for 24 months.                |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Hazard Information</b>  | This product is for research use only, not for diagnostic or therapeutic use.               |
| <b>Dilution Range</b>      | 1:20-1:100 for most applications  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |