

## Product Datasheet

# Affinipure Rabbit Anti-Goat IgG (H+L) antibody (Biotin) (orb633221)

<b>Description</b>	Rabbit polyclonal antibody against Affinipure Rabbit Anti-Goat IgG (H+L) conjugated to Biotin
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Goat
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, ICC, IHC, WB
<b>Preservatives</b>	5mg/ml rAlbumin ( IgG-Free, Protease-Free)
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	5mg/ml BSA ( IgG-Free, Protease-Free)
<b>Storage</b>	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Hazard Information</b>	This product is for research use only, not for diagnostic or therapeutic use.
<b>Dilution Range</b>	1:100 - 1:1000 for enzyme immunohistochemistry on tissue sections 1:4000-1:80000 for ELISA and Western blotting with chromogenic substrates
<b>Expiration Date</b>	12 months from date of receipt.