

Product Datasheet

Affinipure Goat Anti-Rat IgG (H+L) antibody (Biotin) (orb633223)

Description	Goat polyclonal antibody against Affinipure Goat Anti-Rat IgG (H+L) conjugated to Biotin
Species/Host	Goat
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC, WB
Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Concentration	5mg/ml BSA (IgG-Free, Protease-Free)
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:100 - 1:1000 for enzyme immunohistochemistry on tissue sections 1:4000-1:80000 for ELISA and Western blotting with chromogenic substrates
Expiration Date	12 months from date of receipt.