

Product Datasheet

Affinipure Rabbit Anti-Goat IgG (H+L) antibody (Texas Red) (orb633230)

Description	Rabbit polyclonal antibody against Affinipure Rabbit Anti-Goat IgG (H+L) conjugated to Texas Red
Species/Host	Rabbit
Reactivity	Goat
Conjugation	Texas Red
Tested Applications	FC, IF
Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Concentration	5mg/ml BSA (IgG-Free, Protease-Free)
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:20 - 1:100 for most applications
Expiration Date	12 months from date of receipt.