

Product Datasheet

Affinipure Goat Anti-Rabbit IgG (H+L) antibody (Rhodamine) (orb633241)

Description	Goat polyclonal antibody against Affinipure Goat Anti-Rabbit IgG (H+L) conjugated to Rhodamine
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Rhodamine
Tested Applications	FC, IF
Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Concentration	5mg/ml BSA (IgG-Free, Protease-Free)
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:20-1:100
Expiration Date	12 months from date of receipt.