

Product Datasheet

Affinipure Goat Anti Human IgM (H+L) antibody (FITC) (orb633269)

Description	Rabbit polyclonal antibody against Affinipure Goat Anti Human IgM (H+L) conjugated to FITC
Species/Host	Goat
Reactivity	Human
Conjugation	FITC
Tested Applications	FC, IF
Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Concentration	5mg/ml BSA (IgG-Free, Protease-Free)
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	the minimum working dilution was determined to be 1:64 using humanperipheral blood lymphocytes.
Expiration Date	12 months from date of receipt.