

Product Datasheet

Affinipure Rabbit anti Goat IgG(H+L) antibody (Gold) (orb633274)

Description	Rabbit polyclonal antibody against Affinipure Rabbit anti Goat IgG(H+L) conjugated to Gold
Species/Host	Rabbit
Reactivity	Goat
Conjugation	Gold
Preservatives	1%rAlbumin
Form/Appearance	liquid
Storage	2mM boric acid,1%trehalose, 0.1%NaN3, pH9.0, 4°C for 24 months
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Expiration Date	12 months from date of receipt.